


# Vaginal tightness is real.

## The Milli Vaginal Dilator can help.



**Avoid discomfort  
from reinsertion**



#### Indications For Use

The Milli Vaginal Dilator is a tool intended for controlled dilation of the vagina. It can be used for dilation for an examination or to help relieve the symptoms of vaginismus and related dyspareunia.

#### Contraindications

Women with active pelvic infections, or women who have open wounds in the tissue inside or surrounding the vagina should not use vaginal dilators. Pregnant women should refer to their physician for guidance prior to use.

#### Safety Statement

For full safety information, refer to [hellomilli.com](https://hellomilli.com)

#### IFU

For full safety information, including warnings and precautions, refer to [hellomilli.com](https://hellomilli.com)

Learn more:



[hellomilli.com](https://hellomilli.com)

✓ Free Shipping

✓ 1 year warranty

✓ FSA/HSA eligible